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- [1] Liebner S, Cavallaro U & Dejana E. The multiple languages of endothelial cell-to-cell communication. *Arterioscler Thromb Vasc Biol* 26, 1431–1438 (2006).
- [2] Yoon YJ, Kim OY & Gho YS. Extracellular vesicles as emerging intercellular comunicasomes. *BMB Reports* vol. 47 531–539 Preprint at <https://doi.org/10.5483/BMBRep.2014.47.10.164> (2014).
- [3] Tetta C, Ghigo E, Silengo L, Deregibus MC & Camussi G. Extracellular vesicles as an emerging mechanism of cell-to-cell communication. *Endocrine* 44, 11–19 (2013).
- [4] Iraci N, Leonardi T, Gessler F, Vega B & Pluchino S. Focus on extracellular vesicles: Physiological role and signalling properties of extracellular membrane vesicles. *Int J Mol Sci* 17, (2016).
- [5] Yuana Y, Sturk A. & Nieuwland R. Extracellular vesicles in physiological and pathological conditions. *Blood Rev* 27, 31–39 (2013).
- [6] Robbins P D & Morelli AE Regulation of immune responses by extracellular vesicles. *Nat Rev Immunol* 14, 195–208 (2014).
- [7] Raposo G & Stoorvogel W. Extracellular vesicles: Exosomes, microvesicles, and friends. *Journal of Cell Biology* vol. 200 373–383 Preprint at <https://doi.org/10.1083/jcb.201211138> (2013).
- [8] Tamkovich SN., Tutanov OS & Laktionov PP. Exosomes: Generation, structure, transport, biological activity, and diagnostic application. *Biochemistry (Moscow) Supplement Series A: Membrane and Cell Biology* 2016 10:3 10, 163–173 (2016).
- [9] Witwer KW et al. Updating MISEV: Evolving the minimal requirements for studies of extracellular vesicles. *J Extracell Vesicles* 10, (2021).
- [10] Stahl PD & Raposo G. Extracellular Vesicles: Exosomes and Microvesicles, Integrators of Homeostasis. *Physiology (Bethesda)* 34, 169–177 (2019).
- [11] Akers JC, Gonda D, Kim R, Carter BS & Chen CC. Biogenesis of extracellular vesicles (EV): exosomes, microvesicles, retrovirus-like vesicles, and apoptotic bodies. *J Neurooncol* 113, 1 (2013).
- [12] Pan BT & Johnstone RM. Fate of the transferrin receptor during maturation of sheep reticulocytes in vitro: selective externalization of the receptor. *Cell* 33, 967–978 (1983).
- [13] Qiu P, Zhou J, Zhang J, Dong Y & Liu Y. Exosome: The Regulator of the Immune System in Sepsis. *Front Pharmacol* 12, 880 (2021).
- [14] Liu W. et al. Role of Exosomes in Central Nervous System Diseases. *Front Mol Neurosci* 12, 240 (2019).
- [15] Liu J et al. The biology, function, and applications of exosomes in cancer. *Acta Pharm Sin B* 11, 2783–2797 (2021).
- [16] Dai J et al. Exosomes: key players in cancer and potential therapeutic strategy. *Signal Transduction and Targeted Therapy* 2020 5:1 5, 1–10 (2020).
- [17] Hoen EN, Cremer T, Gallo RC & Margolis LB Extracellular vesicles and viruses: Are they close relatives? *Proc Natl Acad Sci U S A* 113, 9155–9161 (2016).
- [18] Lin Y, Anderson JD, Rahnema LMA, Gu Sv & Knowlton AA. Exosomes in disease and regeneration: biological functions, diagnostics, and beneficial effects. *Am J Physiol Heart Circ Physiol* 319, H1162–H1180 (2020).
- [19] Huda MN et al. Potential Use of Exosomes as Diagnostic Biomarkers and in Targeted Drug Delivery: Progress in Clinical and Preclinical Applications. *ACS Biomater Sci Eng* 7, 2106–2149 (2021).
- [20] Xu L, Wu LF & Deng FY. Exosome: An Emerging Source of Biomarkers for Human Diseases. *Curr Mol Med* 19, 387–394 (2019).
- [21] Muthu S, Bapat A, Jain R, Jeyaraman N & Jeyaraman M. Exosomal therapy—a new frontier in regenerative medicine. *Stem Cell Investig* 8, (2021).
- [22] Tejerina S, Dufey V, Hoet JF, Tacheny A & De Longueville F. The DASbox® Mini Bioreactor System as a Tool for Process Development And Stem-Cell Derived Exosome Production in Standardized Culture Conditions. (2021).
- [23] Escobar Ivirico JL, Suttle A & Sha M. SciVario® twin Self-Scale-up of CHO Culture-Based Antibody Production from BioBLU® 3c to BioBLU® 50c Single-Use Bioreactor. www.eppendorf.group/sci-vario (2022).

- [24] Escobar Ivirico JL & Sha M. Stem Cell Exosome Production on the SciVario® twin, a Flexible Controller for Your Bioprocess Needs. (2020).
- [25] Ludwig N, Whiteside TL & Reichert TE. Challenges in Exosome Isolation and Analysis in Health and Disease. *Int J Mol Sci* 20, (2019).
- [26] Sidhom K, Obi PO & Saleem AA. Review of Exosomal Isolation Methods: Is Size Exclusion Chromatography the Best Option? *Int J Mol Sci* 21, 1–19 (2020).
- [27] Tauro BJ et al. Comparison of ultracentrifugation, density gradient separation, and immunoaffinity capture methods for isolating human colon cancer cell line LIM1863-derived exosomes. *Methods* 56, 293–304 (2012).
- [28] Szatanek R, Baran J, Siedlar M & Baj-Krzyworzeka M. Isolation of extracellular vesicles: Determining the correct approach (review). *International Journal of Molecular Medicine* vol. 36 11–17 Preprint at <https://doi.org/10.3892/ijmm.2015.2194> (2015).
- [29] Livshits MA et al. Isolation of exosomes by differential centrifugation: Theoretical analysis of a commonly used protocol. *Sci Rep* 5, (2015).
- [30] Chen J et al. Review on Strategies and Technologies for Exosome Isolation and Purification. *Frontiers in Bioengineering and Biotechnology* vol. 9 Preprint at <https://doi.org/10.3389/fbioe.2021.811971> (2022).
- [31] Kotmakçı M & Akbaba GE. Exosome Isolation: Is There an Optimal Method with Regard to Diagnosis or Treatment? in *Novel Implications of Exosomes in Diagnosis and Treatment of Cancer and Infectious Diseases* (InTech, 2017). doi:10.5772/intechopen.69407.
- [32] Gurunathan S, Kang MH, Jeyaraj M, Qasim M & Kim JH. Review of the Isolation, Characterization, Biological Function, and Multifarious Therapeutic Approaches of Exosomes. *Cells* 8, (2019).
- [33] Momen-Heravi, F et al. Current methods for the isolation of extracellular vesicles. *Biological Chemistry* vol. 394 1253–1262 Preprint at <https://doi.org/10.1515/hsz-2013-0141> (2013).
- [34] Patel GK et al. Comparative analysis of exosome isolation methods using culture supernatant for optimum yield, purity and downstream applications. *Scientific Reports* 2019 9:1 9, 1–10 (2019).
- [35] Gupta S et al. An improvised one-step sucrose cushion ultracentrifugation method for exosome isolation from culture supernatants of mesenchymal stem cells. *Stem Cell Res Ther* 9, 1–11 (2018).
- [36] Lane RE, Korbie D, Anderson W, Vaidyanathan R & Trau M. Analysis of exosome purification methods using a model liposome system and tunable-resistive pulse sensing. *Sci Rep* 5, (2015).

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www.eppendorf.com

Beckman Coulter® and Vi-CELL® are registered trademarks of Beckman Coulter, Inc., USA. BioCoat®, Synthemax® and Corning® are registered trademarks of Corning Inc., USA. RoosterBio® is a registered trademark of RoosterBio Inc., USA. Zetasizer® is a registered trademark of Malvern Panalytical Ltd., UK. Thermo Fisher Scientific® is a registered trademark of Thermo Fisher Scientific Inc., USA. Eppendorf®, the Eppendorf Brand Design, BioBLU®, and ThermoMixer® are registered trademarks of Eppendorf SE, Germany. DASbox® and DASware® is a registered trademark of DASGIP Information and Process Technology GmbH, Germany. All rights reserved, including graphics and images. Copyright © 2024 by Eppendorf SE.

Eppendorf SE reserves the right to modify its products and services at any time. This application note is subject to change without notice. Although prepared to ensure accuracy, Eppendorf SE assumes no liability for errors, or for any damages resulting from the application or use of this information. Viewing the application note alone cannot as such provide for or replace reading and respecting the current version of the operating manual.