


# アングルロータ Angle Rotor T15A45

最高回転速度: 15,000rpm  
 Max. Speed  
 最高遠心加速度: 22,100×g  
 Max. R.C.F.  
 全容量: 60ml=5ml×12  
 Nominal Rotor Capacity



## 適用アダプタ・チューブ・ボトル一覧表 Summary of Adapters, Tubes, Bottles

| アダプタ Adapter<br>品名・Part No.・寸法・材質<br>Name・Part No・Size・Material | チューブ Tube/ボトル Bottle<br>品名・Part No.・寸法・実容量<br>Name・Part No・Size・Actual capacity | 最高回転速度<br>Max. Speed<br>(遠心加速度)<br>(Max. R.C.F)<br>(Rmax) |
|---|---|---|
|   | 5.0ml サンプルングチューブ<br>Sampling tube<br>(市販品) (Marketed) ※2<br>(D)18X(L)60 4.6ml   | ※1<br>15,000rpm<br>(22,100xg)<br>(Rmax:8.80cm)            |
|   | 単位 Unit (mm)  |   |

- ・ ※1: 市販品のチューブ、ボトル等は、メーカー指定の許容回転速度または許容遠心加速度以下で使用してください。許容値を超えて使用するとチューブ、ボトル等が破損する場合があります。
- ・ ※1: When using tubes, bottles, etc. on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tube, bottle, etc. whose allowable speed or RCF is not specified. Otherwise the tubes, bottles, etc. may be broken during operation.
- ・ ※2: 本ロータは、下記に記載したチューブ専用のロータになります。他のチューブを使用しないでください。  
 BIO-BIK 品番 : ST-500
- ・ ※2: Only the tubes given in the following are applicable to this rotor. Do not use tubes other than the following with this rotor.  
 BIO-BIK Catalog No. ST-500