





アングルロータ Angle Rotor T11A31

最高回転速度: 11,000rpm
 Max. Speed
 最高遠心加速度: 11,900×g
 Max. R.C.F.
 全容量: 172ml=43ml×4
 Nominal Rotor Capacity



適用アダプタ・チューブ・ボトル一覧表 Summary of Adapters, Tubes, Bottles

| アダプタ Adapter 品名 ・ Part No. ・ 寸法 ・ 材質 Name ・ Part No ・ Size ・ Material | チューブ Tube/ボトル Bottle 品名 ・ Part No. ・ 寸法 ・ 実容量 ・ 材質 Name ・ Part No ・ Size ・ Actual capacity ・ Material | 最高回転速度 Max. Speed (遠心加速度) (Max. R.C.F) (Rmax) |
|--|--|---|
| _____ |  50ml 培養管 Culture Tube (市販品) (Marketed) (D)30X(L)115 43ml ポリプロピレン Polypropylene | ※1 11,000rpm (11,900xg) (Rmax:8.8cm) |
| _____ |  15ml 培養管 Culture Tube (市販品) (Marketed) (D)17X(L)120 13ml ポリプロピレン Polypropylene | ※1 11,000rpm (11,900xg) (Rmax:8.8cm) |
| <div style="text-align: right;">単位 Unit (mm)</div>  | | |

- ・ ※1: 市販品のチューブ、ボトル等は、メーカー指定の許容回転速度または許容遠心加速度以下で使用してください。許容値を超えて使用するとチューブ、ボトル等が破損する場合があります。
- ・ ※1: When using tubes, bottles, etc. on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tube, bottle, etc. whose allowable speed or RCF is not specified. Otherwise the tubes, bottles, etc. may be broken during operation.
- ・ ※2: このロータは、50ml と 15ml の培養管をアダプタなしで使用できますが、同時には使用できませんので、どちらか一方のみで、使用してください。
- ・ ※2: Either 50-ml or 15ml culture tubes can be used in this rotor without the use of adapters. Do not load 50-ml or 15ml culture tubes together in this rotor for an operation.