

Correspondence table of the rotors and the centrifuges

		Enable		Enable with the updated data		Enable by changing to the updated ROM		Enable by changing to the updated ROM and the converted Door		x		Disable												
Type of Rotor Rotor No. (ml) x (pcs.)	Part No.	Subscription period		Centrifuges					Centrifuges of sales discontinuance															
		on sale	end of sale	CR22H	CR22G	CR21G	CR20G(*1)	CR7	CR26H	CR22F	CR21F	CR20F(*1)	CR22E	CR21E	CR20E(*1)	CR22	CR21	CR20(*1)	CR20B2	CR20B3	SCR20B	SCR18B	CR7B3	CR7B
		2004.08	1999.10	1999.10	1999.10	1999.10	1999.10	1999.10	1985.02	1997.1	1997.1	1997.1	1995.09	1995.09	1995.09	1990.09	1991.10	1992.07	1986.05	1988.09	1981.09	1981.09	1987.03	1992.01

Rotors of sales discontinuance

RR10A 3 0 500ml x 6pcs	90428100	1985.02	1993.03	x	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	7K(4)	10K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	8K(4)	8K(4)	8K(4)	8K(4)	7K(4)	7K(4)		
R10A 3 0 500ml x 6pcs	904308D0	1992.09	2004.10	x	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	7K(4)	10K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	10K(15) 9.5K(4)	8K(4)	8K(4)	8K(4)	8K(4)	7K(4)	7K(4)		
R10A2 3 50/15TC x 12pcs	904337D0	1996.05	2003.12		10K(4)	10K(4)	10K(4)		10K(4)	No.0257 ~ - No.0256 10K(4)	10K(4)	10K(4)	10K(4)	No.0064 ~ - No.0063 10K(4)	No.0062 ~ - No.0061 10K(4)	No.0051 ~ - No.0050 10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	8K(4)	x	x	
R16AF 3 2 250ml x 6pcs	904344D0	1998.10	2002.07	x	13K(4)	13K(4)	13K(4)			No.0266 ~ - No.0265 16K(15) 14K(4)	No.0050 ~ - No.0049 13K(15) 12K(4)	No.0064 ~ - No.0063 13K(15) 12K(4)	No.0062 ~ - No.0061 13K(15) 12K(4)	x	x	x	x	x	x	x	x	x	x	x	
R14AF 1 300ml x 6pcs or 250ml x 6pcs	904332D0	1995.08	1999.10		13K(4)	13K(15) 12K(4)	13K(15) 12K(4)	13K(15) 12K(4)		No.0257 ~ - No.0256 14K(4)	13K(15) 12K(4)	13K(15) 12K(4)	13K(13) 12K(4)	8K(4)	8K(4)	8K(4)	No.0170 ~ - No.0169 13K(15) 12K(4)	No.0688 ~ - No.0687 13K(15) 12K(4)	No.0290 ~ - No.0289 13K(15) 12K(4)	x	x	x	x	x	
R12AF 2 500ml x 4pcs	904333D0	1995.08	2002.07		12K(4)	12K(4)	12K(4)			No.0257 ~ - No.0256 12K(4)	12K(4)	12K(4)	12K(4)	7K(4)	7K(4)	7K(4)	No.0170 ~ - No.0169 12K(15) 11K(4)	No.0688 ~ - No.0687 12K(15) 11K(4)	No.0290 ~ - No.0289 12K(15) 11K(4)	x	x	x	x	x	
R9AF 5 1000ml x 4pcs	904338D0	1997.02	2002.07		9K(4)	9K(4)	9K(4)			x	9K(4)	9K(4)	9K(4)	No.0134 ~ - No.0133 9K(4)	No.0177 ~ - No.0176 9K(4)	No.0099 ~ - No.0098 9K(4)	9K(4)	9K(4)	9K(4)	6K(4)	6K(4)	6K(4)	6K(4)	x	x
R7AF 7 1000ml x 6pcs	90434000	1997.06	1998.08	x	x	x	x	No.0050 ~ - No.0049 7K(15) 6K(4)		x	x	x	x	x	x	x	x	x	x	x	x	x	x	x	
R7AF2 8 1000ml x 6pcs	90434600	1998.08	2002.07	x	x	x	x	No.1004 ~ - No.1003 7K(15) 6K(4)		x	x	x	x	x	x	x	x	x	x	x	x	x	6.6K(15) 5.8K(4)	x	
R15S (RPRS15) 3 5 12ml x 4pcs	90431100	1995.09	1998.05	x	x	x	x			13K(4)	x	x	x	x	x	x	15.0K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(4)	15K(4)	15K(4)	11K(4)	x	x
RPRS15 3 5 12ml x 4pcs	90423900	1989.12	1993.03	x	x	x	x			13K(4)	x	x	x	x	x	x	15.0K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(4)	15K(4)	15K(4)	11K(4)	x	x
R15S2 7 1 12ml x 4pcs	90433600	1995.09	1998.05		15K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(15) 14K(4)		x	15K(15) 14K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(15) 14K(4)	15K(15) 14K(4)	x	x	x	x	x	x	x	x	
RPRS10 1 0 50ml x 4pcs	90423400	1989.02	1992.09	x	10K(4)	10K(4)	10K(4)			10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	10K(4)	8K(4)	x	x
SRR7SA 3 7 250ml x 4pcs	90425100	1983.04	1993.03	x	7K(4)	7K(4)	7K(4)	7K(4)		7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	7K(4)	6K(4)	7K(4)	7K(4)
RR5S2 6 1000ml x 6pcs Bucket aren't attached	90428500	1986.06	1997.03	x	x	x	x	4.2K(4)		x	x	x	x	x	x	x	x	x	x	x	x	x	5K(4)	4.2K(4)	
RPRS4 3 9 50ml x 8pcs	90422800	1986.06	1993.03	x	5K(4)	5K(4)	5K(4)			x	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	5K(4)	x	x
RPRS4-7 7 5ml x 168pcs etc.	90422900	1988.03	1994.03	x	x	x	x	3.2K(4)		x	x	x	x	x	x	x	x	x	x	x	x	x	3.2K(4)	3.2K(4)	
RPRS3-3 4 0 50ml x 16pcs etc.	90423800	1989.01	1993.03	x	3K(4)	3K(4)	3K(4)	3K(4)		x	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)	3K(4)
RPRC18-3 5 4 1000ml	90424600	1981.02	1993.03	x	Batch/ Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)		x	Batch/ Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch x Continuance 18K(15) 16K(4)	Batch/ Continuance No.0100 ~ 18K(15) 16K(4)	Batch/ Continuance 18K(4)	Batch/ Continuance 18K(4)	Batch/ Continuance 18K(4)	Batch/ Continuance 18K(4)	Batch/ Continuance 14K(4)	x	x
RPRC13 5 5 2200ml	90424500	1981.02	1992.09	x	Batch/ Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)		x	Batch/ Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)	Batch x Continuance 13K(4)	Batch/ Continuance 13K(4)	Batch/ Continuance 13K(4)	Batch x Continuance 13K(4)	13K(4)	13K(4)	13K(4)	11K(4)	x	x
RPRC10 5 6 3200ml	90424400	1981.02	1993.09	x	Batch/ Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)		x	Batch/ Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)	Batch x Continuance 10K(4)	Batch/ Continuance 10K(4)	Batch/ Continuance 10K(4)	Batch x Continuance 10K(4)	10K(4)	10K(4)	10K(4)	8K(4)	x	x
R20V2 5 2 5ml x 16pcs	904352D0	1999.10	2004.09		20K(4)	20K(4)	20K(4)			x	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	18K(4)	x	x
R20V 5 1 5ml x 16pcs	90439000	1992.09	1999.10	x	20K(4)	20K(4)	20K(4)			x	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	18K(4)	x	x
RPRV20 5 1 5ml x 16pcs	90423500	1987.01	1993.03	x	20K(4)	20K(4)	20K(4)			x	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	20K(4)	18K(4)	x	x
R11D 6 0 30ml x 8pcs	90431000	1992.09	2004.04	x	11K(4)	11K(4)	11K(4)			x	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	x	x
RPRD11 6 0 30ml x 8pcs	90423700	1991.07	1993.03	x	11K(4)	11K(4)	11K(4)			x	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	11K(4)	x	x

(*1) This centrifuge is sale only in Japan.

(*2) 22K(4) in the table indicates that the max. speed is 22,000min⁻¹ and the control temprature is 4 .

This table shows the correspondence situation to the various centrifuges of the rotor in this time. If you have an unknown point, please ask us freely.